

Biogel® Synthetic Polyisoprene Glove Range



Synthetic Polyisoprene

A sterile, powder-free, synthetic polyisoprene glove that eliminates the possibility of glove-related latex protein sensitization because it is made from a synthetic elastomer. The glove provides levels of fit, feel and comfort comparable to natural rubber latex because synthetic polyisoprene shares a similar molecular structure to natural rubber latex.



Biogel® PI Micro



Biogel® PI UltraTouch®



Biogel® PI UltraTouch® G



Biogel® PI UltraTouch® M



Biogel® PI



Biogel® PI Pro-Fit®



Biogel® PI OrthoPro®



Biogel® PI Indicator® Underglove



Biogel® PI Micro Indicator® Underglove

Description	Enhanced tactile sensitivity	Enhanced sensitivity and softness	Enhanced sensitivity and softness with extra grip	Enhanced sensitivity and softness with extra grip and control	Enhanced softness and medium thickness	Reduces hand and thumb fatigue during surgery	Thicker for increased protection	For use alone or as part of a Puncture Indication System	For use alone or as part of a Puncture Indication System	
Color	Straw Color	Straw Color	Straw Color	Straw Color	Straw Color	Straw Color	Brown Color	Blue for Puncture Indication	Blue for Puncture Indication	
Reorder Number	48555-48590	41155-41190	42155-42190	42655-42690	40855-40890	47960-47990	47660-47690	41655-41690	48955-48990	
Thickness	8.3 mils*	10.6 mils*	10.6 mils*	11.2 mils*	12.2 mils*	10.6 mils*	13.8 mils*	11.0 mils*	8.3 mils*	
Grip	1.5	1.0	1.5	1.5	1.0	1.5	1.5	1.0	1.0	
Finish	Micro-roughened Surface	Micro-roughened Surface	Micro-roughened Surface	Micro-textured Surface	Micro-roughened Surface	Micro-roughened Surface	Micro-roughened Surface	Smooth Finish for Double Donning	Smooth Finish for Double Donning	
Former	Curved Former	Curved Former	Curved Former	Straight Former	Curved Former	Biomechanically Optimized Former	Biomechanically Optimized Former	Curved Former	Curved Former	
Cuff	Anti-Slip, Beaded Cuff	Anti-Slip, Beaded Cuff	Anti-Slip, Beaded Cuff	Beaded Cuff	Anti-Slip, Beaded Cuff	Anti-Slip, Beaded Cuff	Extended, Beaded Cuff	Beaded Cuff	Beaded Cuff	
Sizes	5.5-9.0	5.5-9.0	5.5-9.0	5.5-9.0	5.5-9.0	6.0-9.0	6.0-9.0	5.5-9.0	5.5-9.0	
Quantities	50 Pairs per box/4 boxes per Case 40 pairs per box for size 9	50 Pairs per box/4 boxes per Case 40 pairs per box for size 9	50 Pairs per box/4 boxes per Case 40 pairs per box for size 9	50 Pairs per box/4 boxes per Case 40 pairs per box for size 9	50 Pairs per box/4 boxes per Case 40 pairs per box for size 9	50 Pairs per box/4 boxes per Case 40 pairs per box for size 9	50 Pairs per box/4 boxes per Case 40 pairs per box for size 9	40 Pairs per box/4 boxes per Case	50 Pairs per box/4 boxes per Case 40 pairs per box for size 9	50 Pairs per box/4 boxes per Case 40 pairs per box for size 9
Ideal Use	<i>For use in all surgeries where tactile sensitivity is desired.</i>	<i>For use in all surgeries.</i>	<i>For use in microvascular, cardiothoracic, ENT, plastic, neurosurgical, gynecological and general surgical procedures.</i>	<i>For use in microvascular, cardiothoracic, ENT, plastic, neurosurgical, gynecological and general surgical procedures.</i>	<i>For use in all surgeries.</i>	<i>For use in all surgeries. It can be used alone or as an overglove with Biogel® PI OrthoPro®.</i>	<i>For use in all orthopaedic surgeries, where robust barrier protection is necessary.</i>	<i>For use in all surgeries. Specially designed for ease in double-gloving.</i>	<i>For use in all surgeries. Specially designed for ease in double-gloving and enhanced tactile sensitivity.</i>	

*Finger Thickness Single Wall, typically



The Mölnlycke® Health Care, Biogel®, PI UltraTouch®, Pro-Fit®, OrthoPro®, PI Indicator®, Indicator®, Eclipse®, Super-Sensitive®, Sensor®, Optifit®, NeoDerm® and Skinsense® trademarks, names and logo types are registered globally to one or more of the Mölnlycke® Health Care Group of Companies. Distributed by Mölnlycke® Health Care US, LLC, Norcross, Georgia 30092. © 2016 Mölnlycke® Health Care AB. All rights reserved. 1-800-882-4582. MHC-2015-3847 Lit# 1-0422



Biogel® Latex Glove Range



Natural Rubber Latex

A sterile, powder-free, natural rubber latex glove. All Biogel® latex gloves have exceptionally low extractable protein levels and low endotoxin levels (<0.5 EU/mL).



Biogel® Surgeons



Biogel® Eclipse®



Biogel® Super-Sensitive®



Biogel® M



Biogel® Sensor®



Biogel® Optifit® Orthopaedic



Biogel® Indicator® Underglove



Biogel® Diagnostic



Biogel® NeoDerm®



Biogel® Skinsense®



Biogel® Skinsense® Indicator® Underglove

Description	Provides excellent barrier protection	30% softer with enhanced sensitivity	Enhanced tactile sensitivity with flexibility and manual dexterity	Specially treated surface for extra grip and control	Enhanced tactile sensitivity with a textured surface for a positive grip	10% thicker to provide durability and extra protection	For use alone or as part of a Puncture Indication System	Non-Sterile glove with low dermatitis potential	Low dermatitis potential	Provides extra protection due to longer chemical breakthrough time	For use alone or as part of a Puncture Indication System
Color	Straw Color	Straw Color	Straw Color	Straw Color	Straw Color	Straw Color	Green for Puncture Indication	Straw Color	Khaki Color	Straw Color	Blue for Puncture Indication
Reorder Number	30455-30490	75255-75290	82555-82590	30555-30590	30655-30690	31060-31090	31255-31290	30355-30390	42955-42990	31455-31490	40655-40690
Thickness	10.6 mils*	9.4 mils*	8.5 mils*	11.0 mils*	9.1 mils*	11.8 mils*	8.5 mils*	11.0 mils*	8.7 mils*	8.3 mils*	8.3 mils*
Grip	1.0	2.0	2.0	2.0	2.5	3.0	1.0	1.5	1.5	1.5	1.0
Finish	Micro-roughened Surface	Micro-roughened Surface	Micro-roughened Surface	Micro-textured Surface	Micro-textured Surface	Micro-roughened Surface	Smooth Finish for Double Donning	Micro-textured Surface	Micro-roughened Surface	Micro-roughened Surface	Smooth Finish for Double Donning
Former	Curved Former	Curved Former	Curved Former	Straight Former	Straight Former	Curved Former	Curved Former	Straight Former	Curved Former	Curved Former	Curved Former
Cuff	Beaded Cuff	Anti-Slip, Beaded Cuff	Beaded Cuff	Beaded Cuff	Beaded Cuff	Optifit Beaded Cuff	Beaded Cuff	Beaded Cuff	Beaded Cuff	Beaded Cuff	Beaded Cuff
Sizes	5.5-9.0	5.5-9.0	5.5-9.0	5.5-9.0	5.5-9.0	6.0-9.0	5.5-9.0	5.5-9.0	5.5-9.0	5.5-9.0	5.5-9.0
Quantities	50 Pairs per box/4 boxes per Case 40 pairs per box for size 9	50 Pairs per box/4 boxes per Case 40 pairs per box for size 9	50 Pairs per box/4 boxes per Case 40 pairs per box for size 9	50 Pairs per box/4 boxes per Case 40 pairs per box for size 9	50 Pairs per box/4 boxes per Case 40 pairs per box for size 9	40 Pairs per box/4 boxes per Case	50 Pairs per box/4 boxes per Case 40 pairs per box for size 9	25 Pairs per box/6 boxes per Case	50 Pairs per box/4 boxes per Case 40 pairs per box for size 9	50 Pairs per box/4 boxes per Case 40 pairs per box for size 9	50 Pairs per box/4 boxes per Case 40 pairs per box for size 9
Ideal Use	For use in all surgeries.	For use in all surgeries especially for neurosurgical, plastics and cardiac.	For use in all surgeries especially for neurosurgical, cardiac, transplant and vascular.	For use in ophthalmic, microvascular, cardiothoracic, ENT, plastic, neurosurgical, gynecological.	For use in all surgeries especially for procedures where sensitivity and fingertip control are crucial.	For use in all orthopaedic surgeries where robust barrier protection is necessary.	For use in all surgeries. Specially designed for ease in double-gloving.	For use in all clinical procedures.	For use in all surgeries or where extended protection against chemicals or cytotoxic drugs is required.**	For use in all surgeries or where extended protection against chemicals or cytotoxic drugs is required.**	For use in all surgeries or where extended protection against chemicals or cytotoxic drugs is required.** Specially designed for ease in double-gloving.

*Finger Thickness Single Wall, typically

**Breakthrough testing times are available by request



The Mölnlycke® Health Care, Biogel®, PI UltraTouch®, Pro-Fit®, OrthoPro®, PI Indicator®, Indicator®, Eclipse®, Super-Sensitive®, Sensor®, Optifit®, NeoDerm® and Skinsense® trademarks, names and logo types are registered globally to one or more of the Mölnlycke® Health Care Group of Companies. Distributed by Mölnlycke® Health Care US, LLC, Norcross, Georgia 30092. © 2016 Mölnlycke® Health Care AB. All rights reserved. 1-800-882-4582. MHC-2015-3847 Lit# 1-0422



Biogel® Synthetic Polychloroprene Glove Range



Non-Sterile

Synthetic Polychloroprene

A sterile, powder-free, synthetic polychloroprene glove that reduces the possibility of glove-related latex protein sensitization because it is made from a synthetic elastomer.